



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Armband, Cyan Blue Cloth

**Product Code :**  
LBNY-0012-210008



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Armband, cloth.

Material - 220 gsm 100% washed cotton twill.

The material is color-fast, durable and easy to wash.

Adjustable Velcro closure consisting of two horizontal Velcro strips;  
length 120mm x width 30mm; sewn along the outer side of the armband.

A third vertical Velcro strip sewn on the underside of the armband,  
height 80mm x width 30mm

Color: 100% Cyan Blue

The colour is durable and does not bleach easily in strong sunlight.

There is no leaching of colour when the armband comes into contact with water.

Size: length 420mm x height 80mm.

**Technical Specification :**

Armband, Cyan Blue Cloth



**Laboratory Instrument India**

