



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Beads,wood,col,/BOX-50

**Product Code :**  
LBNY-0013-210002



**Description :**

Beads, wood, not less than 35mm in diameter(+/-10%), with laces, box of 50 pieces, suitable for ages 3 to 6.

**Technical Specification :**

Beads, wooden; not less than 35mm in diameter(+/-10%), in assorted colours (for example: red, green, blue, yellow, orange etc.) and shapes (for example: spheres, cubes, cylinders, ovals, barrel, etc).

The material and finish of the beads and laces should be smooth, secure and non-toxic.

Ten laces for stringing to be provided, of approximately 50cm in length, securely tipped or heat-welded at each end

The design, manufacture and composition of the beads and laces to be suitable for the safe and continued use of children aged 3-6 years.

The laces can be of one or varied colours and must easily fit through the corresponding hole in the shape.

The beads and laces are to be supplied in a sturdy plastic container with a secure lid.

Each bead should have a hole of approximately 6 mm diameter, to enable easy threading of the laces by children.



**Laboratory Instrument India**

