




Email : sales@laboratoryinstrumentindia.com

Product Name :
Professional Trainer Set

Product Code :
LBNY-0002-620001



Description :

Professional Trainer Set, Professional Trainer Set Equipments, Professional Trainer Set Tools, Professional Trainer Set Tool Kits, Professional Trainer Set Manufacturers, Professional Trainer Set Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia

Technical Specification :

Used for the logical linking and galvanic separation of line and domain

In coupler mode the component is wired in parallel on the bus and creates a new line (or zone) to separate the main line and the subline

In repeater mode, the coupler reacts as extender to component line

DUAL CONTROL SWITCH:

Keys can be configured as pairs of (Up and Down) and on/off switch

Individual function programming

Indication LED

Output Switch Actuator, 8 ch, 20A

Power supply, 30V, 38.4Watt, 640mA
Dimmer Actuator, 2ch, N582
Shutter Actuator, 2ch,
USB Interface, dinrail mount
Push Button, 8 fold



Laboratory Instrument India