



Email : sales@laboratoryinstrumentindia.com

Product Name :

Globe,inflatable,diam.42cm,w/o stand

Product Code :

LBNY-0010-210022



Description :

Material: PVC, glossy finish; Approximate thickness 0.19mm; able to survive rough treatment by the children. Inflating valve made of hard plastic, with two plastic loops for inserting rope/thread for hanging. Topographical globe, inflatable (without the use of a pump): showing continents, oceans, rivers, mountains, deserts, forests, et cetera, but without emphasis on political boundaries (printed in dotted lines).

Technical Specification :

Dimensions of globe:

Diameter: 42cm (inflated)

Circumference: approximately 132cm



Laboratory Instrument India