



Email : sales@laboratoryinstrumentindia.com

Product Name :
Compass, plastic, for Blackboard, 45cm

Product Code :
LBNY-0010-2100023



Description :

Compass for blackboard, 40cm long, hard white plastic, with sturdy metal butterfly screw for adjusting tightness by hand. Interchangeable suction foot and compass point in soft plastic. Holder with convenient mechanism that is adjustable for round and square chalk.

Technical Specification :

Weight/Volume/Dimensions:
Estimated volume: 0.82cdm
Estimated weight: 0.227kg



Laboratory Instrument India