




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Clock,teaching,wood,ARABIC

**Product Code :**  
LBNY-0010-210024



**Description :**

Clock; teaching, wood or plastic, with desk-top stand. Minimum diameter of clock 220-240mm; tough construction for classroom handling, easy-to-see numerals, large and small hands.

**Technical Specification :**

Weight/Volume/Dimensions:  
Estimated volume: 5.22 cdm  
Estimated weight: 0.4 kg



**Laboratory Instrument India**