




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Slate, plastic, A4/BOX-20

**Product Code :**  
LBNY-0010-2100027



**Description :**

Suitable for use with HB pencil and eraser (made of natural or synthetic rubber); is washable in water; and sufficiently rigid for the student to be able to write with the slate.

Material: Colour white, matt finish (both sides), high impact polystyrene (HIPS) sheeting, thickness 2mm. resting on the knee.

Slate (individual writing board for students), A4 size (210 x 297mm), white high impact polystyrene sheeting, thickness 2mm, with rounded corners. Silk screen printed on both sides in cyan blue.

Design and rulings: Silk screen printed in cyan blue on both sides.

Reverse side: 8mm lines.

Front side: 1cm squares.

Border: An 8mm border in cyan blue on all sides.

**Technical Specification :**

Weight:

Estimated weight: 2.8kg



**Laboratory Instrument India**

