



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Brush,paint,for chalkboard,50-60mm

**Product Code :**  
LBNY-0010-2100029



**Description :**

The Paint Brush has an excellent application and provides a smooth finish

Brush, paint, for use with two-component paint

Length of bristles: minimum 49mm, maximum 73mm

Width of paintbrush: minimum 50mm, maximum 60mm.

Thickness of brush head: minimum 14mm, maximum 18mm

Synthetic filament preferred.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.329cdm

Estimated weight: 0.1kg



**Laboratory Instrument India**