




Email : sales@laboratoryinstrumentindia.com

Product Name :

Blackboard, double-sided, 1160 x 760 mm

Product Code :

LBNY-0010-2100030



Description :

Chalkboard panel without frame, 1160 x 760 mm, highly pressed fibre board insulite type, smooth surface on both sides, 5 mm (+/-0.5 mm) thick, with rounded corners, finished in black special blackboard paint on both sides and all four edges; granularity suitable for white slate pencil and chalk, UV dried (environmental friendly), for easy cleaning with dry or wet cloth. Two 10mm diameter eyelets (holes) on the 1160 mm side to permit hanging and use of both sides.

Technical Specification :

Weight and Volume:

Estimated volume: 5.405 cdm

Estimated weight: 5.25 kg



Laboratory Instrument India