



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tabards,red/blue,large(age 13+)/PAC-20

**Product Code :**

LBNY-0014-210001



**Description :**

Red and Blue one size fits all coloured tabards to distinguish a sports team, pack of 20 tabards.

**Technical Specification :**

The Size is length 60cm x width 38cm and 10cm velcro on each side. 10 red tabards and 10 blue tabards for a total of 20 tabards. Coloured one size fits all tabards to distinguish a team: made of red/blue, washable, breathable nylon mesh for maximum ventilation, with 10cm adjustable Velcro closures on both sides.



**Laboratory Instrument India**