




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tape, measure, 5 m, retractable

**Product Code :**

LBNY-0014-210002



**Description :**

Tape measure, retractable, for measuring sports grounds; length 5m (five metres), width 13mm; made of thin, plastic-coated metal, calibrated on one side with 1mm and 1cm lengths marked in black and 10cm lengths marked in red, for ease of reading. Stored wound-up in plastic case (assorted colours), 7 x 7cm; with locking mechanism and metal clip, for attachment to referee's clothing.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.13cdm

Estimated dimensions: 7 x 7 x 2.5cm

Estimated weight: 131g



**Laboratory Instrument India**