



Email : sales@laboratoryinstrumentindia.com

Product Name :
Beach Ball, inflatable, diameter 42cm

Product Code :
LBNY-0014-210007



Description :

Inflatable beach ball without the use of a pump with a child friendly design.

Dimensions of beach ball:

Circumference: approx. 132cm

Diameter: approx.42cm (inflated)

Material:

Supplied with:

The repair kit consists of an adhesive vinyl patch measuring 4.5 x 5.5cm.

A transparent vinyl plastic repair kit, with instructions in English for repairing a puncture.

Technical Specification :

with bell inside, diameter 42cm.



Laboratory Instrument India