



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Volleyball, professional model

**Product Code :**  
LBNY-0014-210009



**Description :**

The volleyball is delivered deflated and supplied with two non-return inflating needles. The model is suitable for all types of playground surfaces.

Volleyball, professional/official size, circumference 65-67cm, weight 250-280g. Outer casing is made of synthetic leather, 18 panels; non-return valve; its colour may be a uniform light colour or a combination of colours.

**Technical Specification :**

Weight and Volume:  
Estimated volume: 3.068cdm  
Estimated weight: 319g



**Laboratory Instrument India**