



Email : sales@laboratoryinstrumentindia.com

Product Name :
Volleyball, professional model

Product Code :
LBNY-0014-210009



Description :

The volleyball is delivered deflated and supplied with two non-return inflating needles. The model is suitable for all types of playground surfaces.

Volleyball, professional/official size, circumference 65-67cm, weight 250-280g. Outer casing is made of synthetic leather, 18 panels; non-return valve; its colour may be a uniform light colour or a combination of colours.

Technical Specification :

Weight and Volume:
Estimated volume: 3.068cdm
Estimated weight: 319g



Laboratory Instrument India