



Email : sales@laboratoryinstrumentindia.com

Product Name :
Skipping rope,polyester,3m/PAC-10

Product Code :
LBNY-0014-2100010



Description :

Skipping rope for individual skipping, length 3 metres, diameter 1 centimeter.
Made of soft polyester threads, divided into eight strands, tightly woven, no handles (heat sealed at both ends).
Pack includes ropes of assorted colours.

Technical Specification :

Weight/Volume/Dimensions:
Estimated weight: 0.925kg



Laboratory Instrument India