



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Banner, 150 x 225cm

**Product Code :**  
LBNY-0013-210001



**Description :**

Material: Knitted polyester, at least 110 gsm; 4 eyelets (grommets), one in each corner.  
Size: 150 x 225 cm.

**Technical Specification :**

Weight/Volume/Dimensions:  
Estimated volume: 0.317 cbm  
Estimated weight: 0.77 kg.



**Laboratory Instrument India**