



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Banner,100 x 150cm

**Product Code :**  
LBNY-0013-210002



**Description :**

Material: Knitted polyester, at least 110 gsm; 4 eyelets (grommets), one in each corner.  
Size: 100 x 150 cm.

**Technical Specification :**

Weight/Volume/Dimensions:  
Estimated volume: 0.094 cbm  
Estimated weight: 0.178 kg



**Laboratory Instrument India**