



Email : sales@laboratoryinstrumentindia.com

Product Name :
Flag,outdoors,150 x 225cm

Product Code :
LBNY-0013-210004



Description :

Material: Knitted polyester, at least 110 gsm.

On one side of the flag there is a canvas heading with a rope sewn in (as reinforcement), connecting 2 loops for a flag pole.

Size: 150 x 225 cm.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.207 cbm

Estimated weight: 0.491 kg



Laboratory Instrument India