



Email : sales@laboratoryinstrumentindia.com

Product Name :
First Aid bag,blue,410x280x170mm

Product Code :
LBNY-0013-220001



Description :

First aid bag rectangular, one zippered compartment, padded and lined on top and bottom panels of the compartment with one fixed handle and a detachable shoulder strap.

Shell fabrics: 100% Polyester, 600 x 300 denier, tightly woven and spun to resist puncturing, triple PVC/rubber blend coating. The material measures 340 GSM as a minimum when coated. The product does not permeate by water in the following situations: prolonged water exposure, heavy top-down watering, forceful water pressure applications and submersion.

Lateral fixations and length adjuster: Cyan blue PE plastic.

Cyan piping tape for seams: 100% polyester. Back and front seams of bag reinforced with piping.

Lining fabrics: Cyan taffeta, 100% polyester, 190denier.

Padding: Cyan blue 4mm molitan/sponge rubber, sponge thickness on the sides and top and bottom 10mm. The bag must be padded on all 4 sides and top and bottom.

Zips: Zinc alloy. Zip pull approximately 30mm. 2 Zip pulls Strong zipper with two fasteners sewn on top of the three side panels of the compartment.

Shoulder strap: Cyan blue: 100% polyester. 1200mm x 40mm shoulder strap with rubber pad reinforced with X stitching. Shoulder pad, Cyan blue 4mm rubber.

Carrying Handle: 100% polyester. 30mm wide at the top of the bag. Padded Handle reinforced with X stitching. The bag is reinforced along the side panels with a polyester band measuring 80mm.

Technical Specification :

First Aid bag,blue,410x280x170mm



Laboratory Instrument India