



Email : sales@laboratoryinstrumentindia.com

Product Name :

Bag, blue polyester, 360x230x610mm

Product Code :

LBNY-0013-220002



Description :

The product is not permeated by water in the following situations: prolonged water exposure, heavy top-down watering, forceful water pressure applications, or submersion. Shell fabrics: 100% polyester, 1200 x 600 denier, tightly woven and spun to resist puncturing, triple PVC/rubber blend coating material.

Lateral fixations and Length adjuster

Carrying handle

Cyan piping tape for seams

Shoulder strap

Zips: Zinc alloy

Capacity of Bag: Suitable to carrying a weight of 15kg.

Base reinforcement: Supplied with a loose bottom plate of 1 mm PVC board, dimensions 600 x 230mm (tolerance of +/- 50mm).

Technical Specification :

Carrying bag, rectangular; two handles reinforced with a sponge handle-strap which fastens with Velcro, plus a detachable shoulder strap.

Size: L610 x W230 x H360mm.

Laboratory Instrument India

