



Email : sales@laboratoryinstrumentindia.com

Product Name :

Bag,blue nylon,280x410x170mm

Product Code :

LBNY-0013-220003



Description :

Shell fabrics: 100% Polyester, 1200 x 600 denier, tightly woven and spun to resist puncturing, triple PVC/rubber blend coating material must measure 340 GSM when coated.

Lateral fixations and length adjuster: Cyan blue PE plastic.

Lining fabrics: Cyan taffeta: 100% polyester, 190denier.

Mesh lining: Black 100% polyester. 410x150mm with elastic top and stitched to the base of the bag.

Cyan piping tape for seams: 100% polyester. Back and front seams of bag reinforced with piping.

Shoulder strap: Cyan blue: 100% polyester. 1200 mm x 40 mm shoulder strap with rubber pad reinforced with X stitching. Shoulder pad, Cyan blue 4mm rubber.

Zips: Zinc alloy. Zip pull 30mm.

Padding: 100% polyester, molitan/sponge rubber. The bag is padded on all 4 sides except the top.

Carrying Handle: 100% polyester. 30mm wide at the top of the bag. Handle reinforced with X stitching and padded Velcro closure.

Capacity of Bag: Suitable to carrying a weight of 10kg.

Technical Specification :

Bag,blue nylon,280x410x170mm

Laboratory Instrument India

