



Email : sales@laboratoryinstrumentindia.com

Product Name :
Magnetic Stirrer Hot Plate Magnetic Stir Bar

Product Code :
LBNY-0003-830024



Description :

Magnetic Stirrer Hot Plate Magnetic Stir Bar

Technical Specification :

Magnetic stirrer, Hot Plate, magnetic stir bar

Hot Plate, magnetic stir bar, Magnetic stirrer with hot pot (Max temperature: 280°C) with Detachable foot-long support bar with clamp to allow easy, secure placement of a thermometer or pH probe.

Enclosed cast aluminum heating surface to resist chemical spills.

Sturdy turn dials (one for temperature controls and one for stirring speed controls)

Stirring quantity of H₂O 1 L.

Speed: 0-1600 r.p.m

3pcs of Mixed size (C5/B6 * 15 mm, C6/B7 * 20 mm, C6/B8 * 25 mm/) PTFE magnetic stir bar



Laboratory Instrument India