



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bag,school,students,400x270x100mm

Product Code :
LBNY-0013-220004



Description :

Compartment 1:

Inside Pockets:

Two inner pockets; mesh sides with elastic tops, located on the back side of the main compartment.

Large compartment YKK zip closure around three sides of the compartment with a 2 way YKK zip pull. Mesh lining stitched to the base of the bag and a 1mm PVC board reinforcement secured on the base.

Size: H400 x Depth100 x W270 mm.

Compartment 2:

Small pocket H230 x Depth 40 x W280 mm with side gussets, zip closure around top of pocket with one zip pull.

Technical Specification :

Capacity of Bag: Suitable to carrying a weight of 10kg.



Laboratory Instrument India