



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Bag,school,ARABIC,400x270x100mm

**Product Code :**  
LBNY-0013-220005



**Description :**

Compartment 1:

Inside Pockets:

Two inner pockets; mesh sides with elastic tops, located on the back side of the main compartment.

Large compartment zip closure around three sides of the compartment with a 2 way zip pull. Mesh lining stitched to the base of the bag and a 1mm PVC board reinforcement secured on the base.

Size: H400 x Depth100 x W270 mm.

Compartment 2:

Small pocket H230 x Depth 40 x W280 mm with side gussets, zip closure around top of pocket with one zip pull.

**Technical Specification :**

Bag,school,ARABIC,400x270x100mm

Capacity of Bag: Suitable to carrying a weight of 10kg.



**Laboratory Instrument India**