



Email : sales@laboratoryinstrumentindia.com

Product Name :
Rectangular lapel pin,38 mm long

Product Code :
LBNY-0013-220006



Description :

Lapel pin, metal, rectangular, cyan blue
Material: Stainless steel base, with one butterfly clasp on reverse side. Front: silkscreen printed with two colours and covered with a protective layer of epoxy.
Length of logo: 33mm.
Height of Logo: 4.5mm.
Size: 38x17x10 mm, with rounded corners.

Technical Specification :

Rectangular lapel pin,38 mm long



Laboratory Instrument India