



Email : sales@laboratoryinstrumentindia.com

Product Name :
ID card holder-lanyard/BAG-50

Product Code :
LBNY-0013-220007



Description :

Identity card holder (lanyard), polyester, 20mm wide, in 100% Cyan blue with 1 colour imprint (white) with safety clasp and suspension clip with PVC strap, packed in a bag of 50.

Technical Specification :

The lanyard material is polyester, 20mm in width and 90cm along each length. The polyester is hard wearing, easy to wash, comfortable to wear and with the edges securely finished.

The top of the lanyard contains a plastic safety clasp that is easy to open and close and will release in the event that the lanyard becomes lodged, enabling the wearer to remove the lanyard safely. The clasp is secured at all times.

Identity card holder (lanyard), polyester, in 100% Cyan blue with 1 colour imprint with a plastic safety clasp and fitted with a suspension clip with PVC strap.



Laboratory Instrument India