



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Thermometer, clinical, digital, no Li

**Product Code :**

LBNY-0014-220001



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Clinical thermometer used for taking body temperatures measurements from oral, armpit and rectum.

**INTENDED USE**

A hand-held, battery-powered, electronic instrument designed to measure a patient's body temperature. These variations of the electrical characteristics are processed in the electronic circuits and in turn displayed, for a short period, as temperature readings. Thereafter the display will automatically turn off or go into standby mode. This is a reusable device. It comprises of an electronic unit that detects and converts the changes in temperature into variations of some electrical characteristic, e.g., resistance or voltage.

**Technical Specification :**

Digital, electronic version.  
Glass and mercury free.  
Ready-to-use after switch-on within 10 s.  
Measurement time: within 120 seconds.  
Low and high temperature indication.  
Temperature measurement range 32 – 43 °C (minimum guaranteed).  
Accuracy  $\pm 0.1^{\circ}\text{C}$  in the range 35 – 41 °C.  
Graduation  $0.3^{\circ}\text{C}$  or better.  
Low battery indicator.  
Full batteries allow for a minimum of 4,000 measurements.  
Waterproof.  
Display easy to read in all levels of ambient light.

---

Designed to withstand frequent cleaning and disinfection with hospital-grade products.  
Automatic switch-off when not in use.  
Beep audio alert when device is turned on/ready to use or when temperature measurement is complete.  
Battery powered, batteries included in the supply, preferably packed separately.



**Laboratory Instrument India**