



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Stethoscope,binaural,complete

**Product Code :**  
LBNY-0014-220004



www.laboratoryinstrumentindia.com

**Description :**

It comprises a membrane at the listening head connected by a dual tube to the headgear with ear olives that are placed into the users' ears. A mechanical listening device designed for listening to sounds from the heart, lungs, and/or gastrointestinal tract from adult and paediatric patients.

**Technical Specification :**

Comprises a chest piece connected by a double tube to the headgear with earpieces that are placed into the users' ears.

Designed for frequent and easy disassembly and disinfection with hospital-grade products.

Double cup, with two diaphragms for dual-use (adult and paediatric auscultation) chest piece in zinc alloy.

Tube made of PVC and is crack resistant.

Adult diaphragm Ø: 45,5mm; paediatric diaphragm Ø: 31.5mm.

Tube diameter: outer diameter 10mm, inner diameter 4.8mm. Tube length 560mm.

Sensitivity from 3.2dB to 26dB in a range from 50 to 1000Hz for cardiology.

Sensitivity 8.1dB in a range from 600 Hz to 1,500Hz for pneumology.

Tube impervious to outside noises, guaranteeing full transmission of sound, good auditive quality.

Arms: brass-steel with a flexible spring.

Latex-free.

Removable plastic earpieces.

---

# Laboratory Instrument India

