



Email : sales@laboratoryinstrumentindia.com

Product Name :

Scanner,ultrasound,mobile,w/access

Product Code :

LBNY-0014-220005



www.laboratoryinstrumentindia.com

Description :

Portable ultrasonic scanning systems are used in physician offices, radiology departments, cardiology departments, emergency departments, critical care units and by mobile imaging service providers. The main benefit of point of care ultrasounds units is that it allows physicians to quickly determine whether an abnormality is present so they can make appropriate patient management decisions. They can also have important applications in disaster triage and mass casualty incidents.

Technical Specification :

General purpose ultrasound scanner, with obstetrics package and probes.

System includes a trolley.

The system supports B mode and M mode.

The unit supports user programmable exams (exam pre-sets).

The unit includes: one (1) convex array probe, 2.5-5 MHz and one (1) endovaginal probe, 5.5-7.5 MHz.

An obstetrics analysis package is included which provides for foetal weight and anatomy, gestational age and foetal growth calculations.

Depending on the model supplied Pulsed Wave Doppler (back/white), is included in the unit.

The unit comes standard with packages for: abdomen, gynaecology/obstetrics, cardiology and small parts.

Power requirements: 100 - 240 Volts - 50/60 Hz (not necessarily in a single unit).

Laboratory Instrument India

