



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Long Lighter

**Product Code :**  
LBNY-0003-830025



**Description :**

Long Lighter

**Technical Specification :**

Long Lighter

Filled with liquid-gas

Good quality gas lighters with extra long pipe. Get better reach for your commercial or domestic gas stoves lightening.

This extra long reach gas lighter gives you a safe distance while lightening up big size commercial gas stoves.

Length: 270±30 mm



**Laboratory Instrument India**