



Email : sales@laboratoryinstrumentindia.com

Product Name :
Stethoscope, foetal, Pinard

Product Code :
LBNY-0014-220009



www.laboratoryinstrumentindia.com

Description :

It is designed as a hollow tube (trumpet-shaped) to transmit the foetal heart sound through the inner channel by air conduction and is pressed against the outside of the mother's abdomen and the user will place their head and ear to the listening end of the device. A mechanical listening device that is used for listening to foetal heart sounds.

Technical Specification :

Designed for frequent and disinfection with hospital-grade products.
Small conical instrument, of which the top allows to the practitioner to place the ear, while the flat and circular basis is placed on the abdomen of the pregnant woman.
Monoaural.
Made from unbreakable plastic or aluminium.
Length between: 140 - 180 mm.

WEIGHT AND VOLUME (Packaged)

Volume: 0.6 dm³.

Weight: 50 gr.

Laboratory Instrument India

