



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Cuffs, blood pressure, set 5 sizes

**Product Code :**

LBNY-0014-2200013



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Set of 5 different sizes cuffs to be used with an existing blood pressure meter.

**INTENDED USE**

It is wrapped around the upper arm or the leg (thigh) of the patient. This is a reusable device which may be intended for single-patient or multi-patient reuse. A band-like device that has an inflatable bladder in an inelastic sleeve with one or two connecting tubes with connector (typically locking connectors) that can be connected to a mechanism for inflating and deflating the bladder and/or a sphygmomanometer or a patient monitoring device/system that is used to determine a patient's blood pressure.

**Technical Specification :**

Composed of a cuff containing an inflatable bag.

The inflatable bag has a single tube.

The cuff is fitted with Velcro fastening, enabling a tight and secure fit around arms.

The cuff is made of durable material (e.g. nylon), which is non-deformable, and washable at 30°C.

The cuff is reinforced at both sides.

Material tube: rubber or other suitable material, e.g. silicone rubber, crack resistant.

Set of 5 sizes:

Small adult

Adult

Infant

Child

Large adult

---

Length tube between: 50 to 70 cm.



**Laboratory Instrument India**