



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Flashlight,pre-focused

**Product Code :**  
LBNY-0014-2200014



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Medical-grade pen light, pre-focused  
Diagnostic pen light, generally used to check patient's pupils.

**Technical Specification :**

Reusable model.  
Compact design easily fits in breast pocket; designed with pocket clip.  
Equipped with LED for high intensity light.  
Durable, sturdy plastic or metal housing.  
Suitable for single use or rechargeable batteries: AAA size  
Mechanical switch, which should prevent accidental activation.  
Is able to withstand frequent cleaning and disinfecting with hospital grade products.  
Diameter: 12.5 and 15mm.  
Length: 127 and 134mm.



**Laboratory Instrument India**

