



Email : sales@laboratoryinstrumentindia.com

Product Name :
Flashlight,pre-focused

Product Code :
LBNY-0014-2200014



www.laboratoryinstrumentindia.com

Description :

Medical-grade pen light, pre-focused
Diagnostic pen light, generally used to check patient's pupils.

Technical Specification :

Reusable model.
Compact design easily fits in breast pocket; designed with pocket clip.
Equipped with LED for high intensity light.
Durable, sturdy plastic or metal housing.
Suitable for single use or rechargeable batteries: AAA size
Mechanical switch, which should prevent accidental activation.
Is able to withstand frequent cleaning and disinfecting with hospital grade products.
Diameter: 12.5 and 15mm.
Length: 127 and 134mm.



Laboratory Instrument India

