



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Otoscope, set

**Product Code :**  
LBNY-0014-2200016



www.laboratoryinstrumentindia.com

**Description :**

It consists of a main unit (commonly referred to as the head) with a built-in light source intended to illuminate the interior of the ear canal with a detachable cone-shaped tube (speculum) inserted in the ear canal; it may facilitate the application of air pressure for pneumatic otoscopy. An electrically powered, hand-held device designed for examination of the outer ear canal and tympanic membrane (eardrum) by direct viewing through the ear opening. A dedicated power supply handpiece and/or an insufflation bulb may be included with the head.

**Technical Specification :**

Handle with a built-in LED light source in a mounted head.  
Handle contains on/off switch.  
Including rotatable viewing lens with a magnification between 3 times.  
Includes a port to connect an insufflation bulb.  
Uses 2 x type AA Alkaline batteries as a power source.  
Is able to withstand frequent cleaning and disinfecting with hospital grade products.  
Colour temperature: Cool white in the range 3,200-4,000K.  
Including reusable, autoclavable specula.  
Specula are cleanable with alcohol wipes.

---

# Laboratory Instrument India

