



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tongue depressor, wooden, disp/BOX-100

**Product Code :**

LBNY-0014-2200017



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The device may be used in a healthcare facility, doctor's office, or in the home. This is a single-use device. A hand-held device used to displace the tongue to facilitate examination of the surrounding organs/tissues. It is typically a double-ended flat device with rounded working ends made of wood.

**Technical Specification :**

Finely polished smooth splinter free wood.

Single use

Nonsterile

With rounded extremities for depressing the tongue during examination of the patient's mouth and the throat.

Length: 120-160mm

Width: 14-18mm

Thickness: 1.8-2.2mm



**Laboratory Instrument India**