



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tongue depressor, wooden, disp/BOX-100

Product Code :

LBNY-0014-2200017



www.laboratoryinstrumentindia.com

Description :

The device may be used in a healthcare facility, doctor's office, or in the home. This is a single-use device. A hand-held device used to displace the tongue to facilitate examination of the surrounding organs/tissues. It is typically a double-ended flat device with rounded working ends made of wood.

Technical Specification :

Finely polished smooth splinter free wood.

Single use

Nonsterile

With rounded extremities for depressing the tongue during examination of the patient's mouth and the throat.

Length: 120-160mm

Width: 14-18mm

Thickness: 1.8-2.2mm



Laboratory Instrument India