



Email : sales@laboratoryinstrumentindia.com

Product Name :
Laryngoscope,neonate,set

Product Code :
LBNY-0014-2200021



www.laboratoryinstrumentindia.com

Description :

It has a handle containing batteries to power its light (a small built-in light bulb or fibreoptic light) for airway illumination, and a curved or straight blade of various designs and lengths that can be hinged/interchanged or integral. A hand-held device intended to be used by medical professionals to manipulate the tongue, preventing it from obstructing the oropharynx and enabling a clear view of the trachea for the insertion of an endotracheal (ET) tube prior to the delivery of inhalation anaesthesia and/or ventilation.

Laryngoscope set for neonates, including pre-terms. To manipulate the tongue and enable a clear view of the trachea for surgical/mechanical ventilation/intubation procedures.

Technical Specification :

Includes a cylindrical handle, slightly ribbed, single piece, $\varnothing > 15$ mm
The handle allows fitting of interchangeable blades of different sizes.
Handle is made from non-ferrous material and sealed against ingress of liquids.
Equipped with LED technology light source.
The handle contains a robust on/off switch.
Stud contact, fitting various sizes and types of depressors.
Uses fibre optic for transmission of the light to the tip of the blade.
Blades are able to withstand autoclaving.
With a set of two blades made from stainless steel:
Miller type: Straight Nr 1, length 100 mm \pm 5 mm.
Miller type: Straight Nr 0, length 80 mm \pm 5 mm.



Laboratory Instrument India