



Email : sales@laboratoryinstrumentindia.com

Product Name :

Laryngoscope,adult,child,set

Product Code :

LBNY-0014-2200022



www.laboratoryinstrumentindia.com

Description :

It has a handle containing batteries to power its light (a small built-in light bulb or fibreoptic light) for airway illumination, and a curved or straight blade of various designs and lengths that can be hinged/interchanged or integral. A hand-held device intended to be used by medical professionals to manipulate the tongue, preventing it from obstructing the oropharynx and enabling a clear view of the trachea for the insertion of an endotracheal (ET) tube prior to the delivery of inhalation anaesthesia and/or ventilation.

Technical Specification :

Includes a cylindrical handle, slightly ribbed, single piece, $\varnothing > 25$ mm
The handle contains a robust on/off switch.
The handle allows fitting of interchangeable blades of different sizes.
Handle is made from non-ferrous material and sealed against ingress of liquids.
Blades are able to withstand autoclaving.
With a set of four blades made from stainless steel:
MacIntosh type: Curved Nr 4, length $155\text{mm} \pm 5$ mm.
Miller type: Straight Nr 1, length $100\text{mm} \pm 5$ mm
MacIntosh type: Curved Nr 2, length $110\text{mm} \pm 5$ mm.
MacIntosh type: Curved Nr 3, length $135\text{mm} \pm 5$ mm.
Equipped with LED technology light source.
Stud contact, fitting various sizes and types of depressors.
Uses fibre optic for transmission of the light to the tip of the blade.



Laboratory Instrument India