



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Adult, injection arm trainer

**Product Code :**  
LBNY-0015-220002



**Description :**

Realistic flashback and realistic feel of puncture through artificial skin and vein injection sites. Easily replaceable skin and veins withstanding multiple needle sticks. Realistic simulation of intramuscular and subcutaneous injections with accurate landmarks. Life-size adult arm with realistic anatomical features. Venous access sites available in the antecubital fossa, alongside of forearm and on the dorsum of the hand.

Easily cleanable.

Can provide simulation of: tourniquet position, venepuncture, IV therapy, including by bolus, infusion therapy, blood collection with artificial blood and aseptic techniques.

**Technical Specification :**

Adult, injection arm trainer



**Laboratory Instrument India**