



Email : sales@laboratoryinstrumentindia.com

Product Name :
Anvil 100

Product Code :
LBNY-0004-820016



Description :

Anvil 100

Technical Specification :

This large, modern 100 kg anvil made as a result of close cooperation between Polish and Swiss.

It is the best for those who are interested in artistic metalworking.

This single horn anvil has a hardy hole for attaching mandrel tools and three pritchel holes of different diameter to punch holes in the material.

The horn of the Perun anvil is a slim cone in shape with a shifted symmetry axis.

The special shape of the horn allows the most complex but precise work such as bending material at any angle or making the rings with steady, regular diameter.

The face of the Perun anvil is divided into two parts: a convex one, which is used to draw the material or to bend it, and a flat one, which can be useful in other blacksmith works.

It is important that both parts are based on anvil's solid body so the energy from forging is well absorbed.

Face: 10.3cm*14.8cm

Heel: 7.8cm*12.7cm

Horn: 8.5cm*21cm

Height: 21cm

Base: 28,2cm*22.5cm

Hammer: 10*4*27cm



Laboratory Instrument India