



Email : sales@laboratoryinstrumentindia.com

Product Name :
Adult, lumbar puncture trainer

Product Code :
LBNY-0015-220003



Description :

Epidural and spinal injection simulator to perform both lumbar puncture and epidural procedures.
Realistic tactile feedback, epidural space, dura, subarachnoid space with fluid supply system
Reliable palpable anatomical landmarks (iliac crests and spinous processes L3/L4 or L4/L5) for orientation of user to access points
Lower torso segment with multilayer lumbar puncture inserts
Can be positioned in both the upright and lateral decubitus positions for various training scenarios
Realistic tissue response including loss of resistance when entering the epidural space
Needle insertion possible between vertebrae to the epidural space for injection of local anaesthetic or access into the subarachnoid space to perform a 'spinal tap'
Durable, easy to clean and maintain
For training in aseptic techniques at puncture site.

Technical Specification :

Supplied with:
1x base body
1x spinal canal with connector
1x water collection container
1x10 cm tube
1x vertebral bar
1x LOR insert (loss-of-resistance)
1x spacer ring for injection
1x extension tube for the hanging-drop technique

1x skin
2x 100 ml injections with cap
1x CD with instructions
1x carry bag



Laboratory Instrument India