



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Newborn,advanced care trainer

**Product Code :**  
LBNY-0015-220004



**Description :**

Full term newborn size manikin with soft, life-like skin and moulded hair. Transparent chest allows direct observation of ribs, heart and superior vena cava.

Facilitates training of following skills:

Central venous cannulation through cephalic and saphenous vein punctures.

Tube feeding and gastric suction.

Translucent skin allows visualization of underlying veins.

Correct anatomical position of basilic, cephalic, jugular and subclavian veins.

**Technical Specification :**

Supplied with:

One extra vein with attached tubing

One bottle of simulated blood

One right arm skin

One right leg skin

Two nasogastric tubes.



**Laboratory Instrument India**

