



Email : sales@laboratoryinstrumentindia.com

Product Name :
Newborn, injection arm trainer

Product Code :
LBNY-0015-220006



Description :

Life-size newborn arm. With realistic anatomical features such as: medial venous structures; veins on the dorsum of the hand; radial and brachial arteries.
Easily replaceable resealing skin and veins withstanding multiple needle sticks. Realistic “pop” when needle enters vein.
Squeeze bulb to simulate collapsed and normal veins.

Technical Specification :

Supplied with:
One squeeze bulb
One spare skin for arm
One bag for artificial blood with stand
One synthetic blood concentrate
Carrying bag.



Laboratory Instrument India

