



Email : sales@laboratoryinstrumentindia.com

Product Name :
Fetal baby,w/umbilical cord & placenta

Product Code :
LBNY-0015-220008



Description :

Foetal baby, with umbilical cord and placenta.
Full size, flexible term foetus with lifelike features
Umbilical cord and placenta present for training in cord ligation
Realistic placenta for training in placenta expulsion exercises
With smooth, soft skull containing anterior and posterior fontanelles for training in realistic delivery,
Can be used in training of Vacuum deliveries
Should fit into advanced childbirth simulator.

Technical Specification :

Weight and volume:
Estimated Volume: 0.03 m3
Estimated Weight: 2 kg



Laboratory Instrument India