



Email : sales@laboratoryinstrumentindia.com

Product Name :
Labor delivery module

Product Code :
LBNY-0015-220009



Description :

The labour delivery module consists of six selected labour stations to represent conditions of the cervix and vagina prior to labour, during labour, and at birth in a primigravida
It provides training in assessing the position of the fetal head in relation to the maternal pelvis and the level of the ischial spine

Stations illustrated are in this model are:

- STA -3 cervix fully effaced;
- STA 0 foetal head at plane of ischial spine;
- STA +2 cervix nearing full dilation;
- STA -5 prior to onset of labour;
- STA -4 cervix partially effaced;
- STA +5 crowning of foetal head

Stations may be placed in the birth canal of the advanced childbirth simulator
Should fit into Advanced childbirth simulator.

Technical Specification :

Weight and volume:
Estimated Volume: 0.05 m³
Estimated Weight: 6.5kg

Laboratory Instrument India

