



Email : sales@laboratoryinstrumentindia.com

Product Name :
Postpartum suturing trainer

Product Code :
LBNY-0015-2200010



Description :

Postpartum suturing simulator that provides training on suturing episiotomies, labial and perineal tears
Realistic palpation of muscular structures including M. sphincter ani externus, M. transversus perinei superficialis and M. bulbospongiosus
Muscle and skin layers can be sutured separately
Allows training on the correct amount of tension to apply to sutures and different knot-tying methods
Provides realistic training on a variety of surgical repair techniques
Soft skin with lifelike texture to provide a realistic feel for suturing and instrument handling
With the possibility of practicing 3 perineal repair procedures:
Midline episiotomy
Standard mediolateral episiotomy
Mediolateral episiotomy.

Technical Specification :

Weight and volume:
Estimated Volume: 0.007 m³
Estimated Weight: 1kg

Laboratory Instrument India

