



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pipe Vice 900

**Product Code :**  
LBNY-0004-820017



**Description :**

Pipe Vice 900

**Technical Specification :**

Pipe Vice 900  
Quick release mechanism secures pipe in no time.  
Simple bolt-on installation for workbenches and work stands.  
Pipe threader plumbing vise tool  
Size: 1/2" - 4-1/2" (10-90mm)  
Material: stainless steel.



**Laboratory Instrument India**