



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Newborn basic care trainer, dark

**Product Code :**  
LBNY-0015-2200013



**Description :**

Material of the simulator: Thermoplastic polyurethane - TPU (primary material for body skin), Poly vinyl chloride - PVC (primary material in chest and face skin) silicone rubber, Polyoxymethylene - POM, polypropylene - PP, and iron/steel.

Realistically sculpted to resemble and simulate a baby, anatomically accurate and feeling smooth / resilient / bony as relevant and suitable for simulation. Full term anatomically accurate newborn with umbilical cord

Texture of the simulator: soft and smooth, friction free to demonstrate the desired procedures.

**Technical Specification :**

Environmental conditions:

Storage conditions: -10°C to +50°C / 5-95% RH. (Short term for upper temperature limit - 2 weeks)

Operating conditions:

Operating humidity: 5 to 95% RH

Operating temperature: +15°C to +40°C



**Laboratory Instrument India**

