



Email : sales@laboratoryinstrumentindia.com

Product Name :
Newborn basic care trainer, light

Product Code :
LBNY-0015-2200014



Description :

Material of the simulator: Thermoplastic polyurethane - TPU (primary material for body skin), Polyvinylchloride - PVC (primary material in chest and face skin) silicone rubber, Polyoxymethylene - POM, polypropylene - PP, and iron/steel.

Full term anatomically accurate newborn with umbilical cord.

Texture of the simulator: soft and smooth, friction free to demonstrate the desired procedures. Realistically sculpted to resemble and simulate a baby, anatomically accurate and feeling smooth / resilient / bony as relevant and suitable for simulation.

Technical Specification :

Environmental conditions:

Storage conditions: -10°C to +50°C / 5-95% RH. (Short term for upper temperature limit - 2 weeks)

Operating conditions:

Operating humidity: 5 to 95% RH

Operating temperature: +15°C to +40°C



Laboratory Instrument India

